

UNIVERSITY OF CALIFORNIA, MERCED
SCIF CLEANROOM FACILITY
STANDARD OPERATING PROCEDURE (SOP)
Atomic Force Microscopy (AFM) Operation (Solver NEXT or Equivalent)
Location: SE1: 153

1. PURPOSE

To provide a standardized and safe procedure for nanoscale surface imaging and characterization using Atomic Force Microscopy (AFM), ensuring high-quality data, instrument protection, and user safety.

2. SCOPE

Applicable to all trained users performing AFM measurements including contact, semicontact (tapping), and non-contact modes for surface imaging and analysis.

3. RESPONSIBILITIES

- **Users:** Follow SOP and operate within approved parameters
 - **Core Staff:** Maintain system and provide training
 - **Facility:** Ensure environmental stability and compliance
-

4. SYSTEM OVERVIEW

The AFM is a scanning probe microscope that measures surface topography using a sharp probe interacting with the sample surface.

Capabilities:

- Nanoscale topography imaging
 - Surface roughness measurement
 - Phase and material contrast imaging
 - Force-distance measurements
 - Advanced modes (MFM, EFM, etc., depending on configuration)
-

5. REQUIRED SYSTEM CONDITIONS

5.1 Utilities

- AC Power (110/220 V)
- AFM controller and software operational
- Vibration-isolated table
- Stable temperature environment (20–25°C recommended)

5.2 System Conditions

- Probe (cantilever) properly installed
- Sample securely mounted
- Scanner calibrated
- Laser aligned on photodiode
- Minimal vibration and airflow disturbance

6. SAFETY REQUIREMENTS

6.1 Hazards

- Laser radiation exposure
- Fragile probe (tip damage risk)
- Piezo scanner damage
- Electrical components

6.2 PPE

- Cleanroom gloves
- Safety glasses
- Lab coat

6.3 Critical Safety Rules

- Do NOT stare into laser beam
- Do NOT apply excessive force to scanner

- Avoid mechanical shock to system
 - Do NOT touch probe tip
 - Ensure proper grounding
-

7. CLEANROOM PROTOCOL

- Handle samples with tweezers only
 - Ensure samples are clean and dry
 - Avoid dust, fibers, and residues
 - Use clean sample holders
 - Minimize exposure to ambient air
-

8. DETAILED OPERATION PROCEDURE

STEP 1: SYSTEM PRE-CHECK

- Verify:
 - o System powered ON
 - o Software running
 - o Probe installed correctly
 - o Sample mounted securely
-

STEP 2: LASER ALIGNMENT

- Align laser on cantilever
 - Adjust photodiode signal
 - Optimize deflection signal
-

STEP 3: SELECT MODE

- Choose imaging mode:
 - o Contact mode
 - o Semicontact (tapping) mode
 - o Non-contact mode
-

STEP 4: APPROACH SAMPLE

- Use automated or manual approach
 - Slowly bring probe toward sample
 - Monitor deflection signal
-

STEP 5: SET SCANNING PARAMETERS

- Scan size (e.g., 1–100 μm)
 - Scan rate (0.5–2 Hz typical)
 - Feedback gain
 - Setpoint (force/amplitude)
-

STEP 6: PERFORM SCAN

- Start scanning
 - Monitor:
 - o Topography signal
 - o Feedback stability
 - o Noise levels
-

STEP 7: OPTIMIZE IMAGE

- Adjust:
 - o Gain
 - o Setpoint
 - o Scan speed
 - Minimize noise and artifacts
-

STEP 8: SAVE DATA

- Save scan images
 - Export data if needed
 - Record in logbook
-

STEP 9: WITHDRAW PROBE

- Retract probe safely from surface
-

STEP 10: REMOVE SAMPLE

- Remove sample carefully
 - Avoid contact with probe
-

STEP 11: SHUTDOWN

- Close software
 - Leave system in safe standby condition
-

9. TROUBLESHOOTING GUIDE

Issue	Cause	Action
No signal	Laser misalignment	Realign laser
No image	Improper approach	Re-approach
Noisy data	Vibration	Stabilize environment
Tip crash	Aggressive approach	Reduce speed
Poor resolution	Damaged tip	Replace probe

10. CRITICAL DOs & DON'Ts

DO

- Align laser properly
- Use appropriate scan parameters
- Maintain vibration isolation
- Monitor scan continuously

DON'T

- Crash probe into sample
- Touch probe tip

- Operate under unstable conditions
 - Use contaminated samples
-

11. CONTAMINATION CONTROL

- Use clean substrates
 - Avoid outgassing materials
 - Keep sample stage clean
 - Store probes properly
-

12. WASTE HANDLING

- Dispose wipes per cleanroom protocol
 - Follow SCIF waste procedures
-

13. EMERGENCY PROCEDURES

- Probe crash → Stop immediately
 - Laser issue → Shut down system
 - System error → Abort scan
 - Electrical issue → Power OFF
-

14. TRAINING COVERAGE

Users are trained on:

- AFM fundamentals
 - Probe handling and installation
 - Imaging modes
 - Parameter optimization
 - Data analysis
-

15. ACKNOWLEDGMENT & APPROVAL

Director Name: _____

Director Signature: _____

Date: _____

User Name: _____

User Signature: _____

Date: _____