

**UNIVERSITY OF CALIFORNIA, MERCED**

**SCIF CLEANROOM FACILITY**

**STANDARD OPERATING PROCEDURE (SOP)**

**Atomic Layer Deposition (ALD) System Operation (Savannah S100/S200/S300 or Equivalent)**

**Location: SE1: 154, Class 100/1000 Cleanroom**

---

## **1. PURPOSE**

To provide a detailed, standardized, and safe procedure for thin film deposition using Atomic Layer Deposition (ALD), ensuring reproducibility, conformality, and precise thickness control while protecting equipment and users.

---

## **2. SCOPE**

Applicable to all trained users performing thin film deposition using ALD systems for oxides, dielectrics, and nanostructured coatings.

---

## **3. RESPONSIBILITIES**

- Users: Follow SOP and ensure proper recipe parameters
  - Core Staff: Maintain system and provide training
  - Facility: Ensure utilities, safety, and compliance
- 

## **4. SYSTEM OVERVIEW**

Atomic Layer Deposition (ALD) is a cyclic thin film deposition technique based on alternating, self-limiting surface reactions.

Typical ALD cycle:

- Precursor A pulse
- Purge
- Precursor B pulse
- Purge

Capabilities:

- Atomic-scale thickness control ( $\sim\text{\AA}$  per cycle)

- Highly conformal coatings
  - Uniform deposition on high aspect ratio structures
  - Thin films from nm to tens of nm scale
- 

## **5. REQUIRED SYSTEM CONDITIONS**

---

### **5.1 Utilities**

- Nitrogen gas ( $\geq 99.999\%$  purity recommended)
  - Vacuum pump operational ( $\sim 10^{-1}$  to  $10^{-3}$  Torr typical base range)
  - System heaters operational (chamber, lines, precursor)
  - Electrical supply (115 VAC system)
- 

### **5.2 System Conditions**

- Chamber clean and free of contamination
  - Precursors properly installed and sealed
  - No vacuum leaks
  - Temperature controllers functioning
  - Proper gas flow and pressure control
- 

## **6. SAFETY REQUIREMENTS**

---

### **6.1 Hazards**

- High temperature surfaces (150–250°C)
  - Hazardous chemical precursors (some pyrophoric)
  - Vacuum system hazards
  - Pressurized gas systems
  - Heated gas lines and valves
- 

### **6.2 PPE**

- Cleanroom gloves
  - Safety glasses
  - Lab coat
- 

### 6.3 Critical Safety Rules

- Never open chamber below operating temperature (~150°C)
  - Never disconnect precursor cylinders improperly
  - Ensure nitrogen flow before operation
  - Do not operate without vacuum
  - Follow all SDS guidelines for precursors
- 

## 7. DETAILED OPERATION PROCEDURE

---

### STEP 1: SYSTEM PRE-CHECK

- Verify:
    - o Nitrogen supply ON
    - o Precursors available
    - o Chamber clean
    - o No leaks
    - o System powered
- 

### STEP 2: SYSTEM STARTUP

- Click **PUMP**
  - Turn **HEATER ON**
  - Allow chamber to reach ~150°C (~15 minutes)
- 👉 Do NOT open chamber before temperature is reached
- 

### STEP 3: LOAD SAMPLE

- Click **VENT**
- Open chamber lid
- Place substrate (face up) using tweezers

- Close chamber
  - Click **PUMP**
  - Allow thermal stabilization (~5–10 minutes)
- 

#### **STEP 4: SET PROCESS CONDITIONS**

- Set nitrogen flow: ~5–20 sccm
  - Confirm pressure (~100 mTorr typical)
  - Verify temperature (~150°C or as required)
- 

#### **STEP 5: RECIPE SETUP**

- Load existing recipe template
  - Modify:
    - o Number of cycles
    - o Pulse times
    - o Purge times
  - Save recipe under user directory
- 👉 Typical growth rate: ~1 Å per cycle (material dependent)
- 

#### **STEP 6: RUN PROCESS**

- Click **RUN (START)**
  - Confirm execution
  - Monitor:
    - o Pulse signals
    - o Pressure
    - o Cycle progression
- 

#### **STEP 7: PRECURSOR MONITORING**

- Ensure pulses are present and consistent
- 👉 If precursor pulse missing → check valve
- 👉 If water pulse missing → STOP and notify staff

---

### **STEP 8: DEPOSITION**

- System runs automatically through cycles
- Monitor:
  - o Pressure stability
  - o Pulse amplitude
  - o Temperature

---

### **STEP 9: COMPLETE PROCESS**

- Wait for “RUN HAS COMPLETED”
- Record process in logbook

---

### **STEP 10: UNLOAD SAMPLE**

- Click **VENT**
- Open chamber
- Remove sample using tweezers
- Close chamber

---

### **STEP 11: PURGE MANIFOLD (MANDATORY)**

- Close precursor valves
  - Load purge recipe
  - Run purge cycle
- 👉 Repeat for all precursors used

---

### **STEP 12: SHUTDOWN**

- Maintain low nitrogen flow (~5 sccm)
- Leave system in standby
- Do NOT shut down software unless emergency

---

## **8. TROUBLESHOOTING GUIDE**

<b>Issue</b>	<b>Cause</b>	<b>Action</b>
No precursor pulse	Valve closed	Open valve
No water pulse	Empty source	Notify staff
Pressure instability	Leak / contamination	Stop process
Poor film quality	Incorrect parameters	Adjust recipe
Non-uniform film	Temperature/flow issue	Verify conditions

---

## **9. CRITICAL DOs & DON'Ts**

### **DO**

- Use proper purge cycles
- Maintain clean chamber
- Monitor process continuously
- Use approved precursors

### **DON'T**

- Open chamber when cold or mid-process
  - Skip purge steps
  - Use unknown precursors
  - Operate with leaks or unstable pressure
- 

## **10. CONTAMINATION CONTROL**

- Clean chamber regularly
  - Avoid outgassing substrates
  - Use clean handling tools
  - Maintain precursor line cleanliness
- 

## **11. WASTE HANDLING**

- Dispose of wipes and materials properly
- Follow SCIF hazardous waste protocols
- Handle precursor waste per SDS

---

## 12. EMERGENCY PROCEDURES

- Power issue → Shut OFF immediately
  - Vacuum failure → Stop process
  - Gas issue → Close supply
  - Abnormal behavior → Abort run
- 

## 13. TRAINING COVERAGE

Users are trained on:

- ALD fundamentals
  - Recipe programming
  - Precursor handling
  - System operation
  - Safety procedures
- 

## 14. ACKNOWLEDGMENT & APPROVAL

Director Name: \_\_\_\_\_

Director Signature: \_\_\_\_\_

Date: \_\_\_\_\_

User Name: \_\_\_\_\_

User Signature: \_\_\_\_\_

Date: \_\_\_\_\_