

UNIVERSITY OF CALIFORNIA, MERCED

SCIF CLEANROOM FACILITY

STANDARD OPERATING PROCEDURE (SOP)

UV Exposure System Operation (350W Collimated NUV Exposure System or Equivalent)

Location: SE1: 154, Class 100 Cleanroom

1. PURPOSE

To provide a standardized and safe procedure for ultraviolet (UV) exposure processes used in photolithography, ensuring uniform exposure, reproducibility, and user safety.

2. SCOPE

Applicable to all trained users performing UV exposure for photoresist patterning, flood exposure, resist stabilization, and related lithography processes.

3. RESPONSIBILITIES

- **Users:** Follow SOP and exposure parameters
 - **Core Staff:** Maintain system and provide training
 - **Facility:** Ensure safe operation and compliance
-

4. SYSTEM OVERVIEW

The UV exposure system uses a high-intensity mercury arc lamp and optical system to generate a uniform, collimated UV beam for photolithography.

Capabilities:

- Photoresist exposure (350–450 nm NUV range)
- Uniform flood exposure
- Edge bead exposure
- High uniformity ($\pm 5\text{--}6\%$)

👉 The system uses dielectric mirrors and an optical integrator to produce a uniform beam

5. REQUIRED SYSTEM CONDITIONS

5.1 Utilities

- AC Power supply
 - UV lamp power supply operational
 - Cooling/exhaust system active
 - Shutter control system functional
-

5.2 System Conditions

- Lamp properly installed
 - Optical components clean
 - Shutter functional
 - System covers in place
 - Proper ventilation available
-

6. SAFETY REQUIREMENTS

6.1 Hazards

- Intense UV radiation (eye/skin damage)
 - Lamp explosion risk
 - Ozone generation
 - Electrical shock (high voltage)
 - Mercury exposure (lamp)
-

6.2 PPE

- UV-blocking safety glasses
 - Cleanroom gloves
 - Lab coat
 - Face shield (during lamp handling)
-

6.3 Critical Safety Rules

- NEVER look directly at UV beam
- ALWAYS wear UV eye protection
- Do NOT operate without covers in place
- Ensure proper ventilation (ozone removal)
- Disconnect power before maintenance

👉 UV radiation and lamp explosion are major hazards requiring strict precautions

7. CLEANROOM PROTOCOL

- Handle wafers using tweezers
 - Ensure photoresist-coated samples are clean
 - Avoid contamination on exposure stage
 - Keep optics and stage particle-free
-

8. DETAILED OPERATION PROCEDURE

STEP 1: SYSTEM PRE-CHECK

- Verify:
 - o Power supply ON
 - o Cooling/exhaust active
 - o Lamp installed correctly
 - o Covers secured
-

STEP 2: WARM-UP

- Turn ON UV lamp
- Allow warm-up (10–15 minutes)

👉 Lamp must reach stable operating conditions for uniform exposure

STEP 3: SET EXPOSURE PARAMETERS

- Set lamp power (typically 280–380 W for 350W system)
 - Set exposure time
 - Confirm intensity settings
-

STEP 4: LOAD SAMPLE

- Place wafer/sample on exposure stage
 - Align sample properly
 - Ensure flat contact
-

STEP 5: ALIGN AND VERIFY INTENSITY

- Place UV meter (if required)
- Verify uniform intensity
- Adjust lamp Z-position for peak intensity

👉 “Peaking” the lamp maximizes exposure intensity

STEP 6: PERFORM EXPOSURE

- Close system cover
 - Activate shutter
 - Expose for defined time
-

STEP 7: COMPLETE EXPOSURE

- Close shutter
 - Remove sample
 - Proceed to development process
-

STEP 8: SHUTDOWN

- Turn OFF lamp (if not in continuous use)
 - Allow system to cool
 - Turn OFF power supply if required
-

9. TROUBLESHOOTING GUIDE

| Issue | Cause | Action |
|----------------------|-------------------|----------------------|
| Low intensity | Lamp aging | Replace lamp |
| Non-uniform exposure | Misalignment | Re-align lamp |
| Overheating | Poor cooling | Check exhaust |
| Ozone smell | Ventilation issue | Increase airflow |
| No exposure | Shutter failure | Check control system |

10. CRITICAL DOs & DON'Ts

DO

- Allow proper lamp warm-up
- Use UV protection
- Maintain clean optics
- Monitor intensity regularly

DON'T

- Look into UV beam
 - Operate without covers
 - Touch lamp quartz surface
 - Adjust optics unnecessarily
-

11. CONTAMINATION CONTROL

- Keep optics clean (lens, mirrors)
 - Avoid touching optical surfaces
 - Clean with IPA and lint-free wipes
 - Maintain clean exposure stage
-

12. WASTE HANDLING

- Dispose lamps as hazardous waste (mercury-containing)
 - Follow SCIF hazardous waste procedures
-

13. EMERGENCY PROCEDURES

- Lamp explosion → Evacuate area, allow cooling
 - UV exposure → Seek medical attention
 - Electrical issue → Power OFF immediately
 - Ozone buildup → Ventilate area
-

14. TRAINING COVERAGE

Users are trained on:

- UV lithography principles
 - Exposure dose control
 - Lamp handling and safety
 - System calibration and intensity measurement
 - Emergency procedures
-

15. ACKNOWLEDGMENT & APPROVAL

Director Name: _____

Director Signature: _____

Date: _____

User Name: _____

User Signature: _____

Date: _____